


A validated SBSE-GC-MS/MS method for ultra-trace analysis of priority pollutants in surface water in the context of the EU water framework directive with additional focus on suspended matter



Oliver Lerch, Jasmin Zboron
GERSTEL GmbH & Co. KG
Muelheim, Germany
Oliver_Lerch@gerstel.de
September 24th 2019



GERSTEL



~0,4 ng/L im Baldeneysee

Joint Project of **Agilent** and **Gerstel**

GERSTEL



Many thanks to the following Agilent colleagues

Juan-Luis Aybar

John Upton

Tarun Anumol

Ken Brady

Alain Vervaecke

Many thanks to the following GERSTEL colleague

Jasmin Zboron



Stuart Miles, freedigitalphotos.net


GERSTEL




Agilent

Trusted Answers

Background on EU Water Framework Directive (EU-WFD)



Priority Compounds(-Groups) EU-WFD




■ GC, included
 ■ GC, not included
 ■ VOCs
 ■ LC
 ■ Metal


- ▷ **GC compatible, included ~70% of compounds:**
 - Around 100 compounds
32 compound groups & individual compounds from EU-WFD
Country specific relevant compounds for surface water analysis
- ▷ **GC compatible, not included:**
 - C10-C13 Chloroalkanes:
dedicated method needed
 - Tributyltin:
dedicated method needed (e.g. RIC Belgium)
 - Dioxins:
could be included, but not in our lab
- ▷ **VOCs:**
 - Dedicated headspace GC-MS/MS method to be developed
- ▷ **LC compatible:**
 - Diuron, isoproturon, hexabromocyclododecane, PFOS
- ▷ **Metals:**
 - With AAS or ICP-MS




Choice of Analytes from 2013/39/EU




Dieldrin




p,p'-DDT




Chlorpyrifos




PCBs




Alachlor




Atrazine




4-tert-Octylphenol




PBDEs




Benzo[a]pyrene



Diethylhexylphthalate



Dichlorvos



1,2,4 Trichlorobenzene

Challenges

- ▶ Numerous differing compounds
- ▶ Low required LOQs
derived as 0.3x the respective annual average EQS: e.g. 600 ng/L naphthalene, p,p'-DDT 3 ng/L, endosulfan 1.5 ng/L, dichlorvos 0.18 ng/L, benzo[a]pyrene 0.051 ng/L, 0.024 ng/L cypermethrine, 0.06 pg/L heptachlor and heptachlor epoxide
- ▶ “Where those criteria are not met for any matrix, Member States shall ensure that monitoring is carried out using best available techniques not entailing excessive costs”
- ▶ “Concentrations in the whole water sample”
-> Extraction of particle adsorbed compounds required

GERSTEL



Agilent

Trusted Answers

Stir Bar Sorptive Extraction (SBSE) GERSTEL Twister (PDMS)


Twister Extraction




1 cm Twister



Add conditioned Twister to sample vial



Parallel extraction – partitioning of the analytes between water and PDMS phase



Removal of Twister after extraction



Removal of residual sample and particles by dipping into clean water



Transfer to clean TDU desorption tube



Drying in N₂ stream, alternatively by clean and lint free tissue www.gerstel.com



Ready to put into autosampler tray

Thermal Desorption Unit
TDU

GC carrier gas

100 % analyte transfer
onto column is possible

Cooled Injection System as
„Cryotrap“
CIS

GERSTEL



Agilent

Trusted Answers

Optimization of Sample Preparation Method

Method Optimization – TDU/CIS Parameters

GERSTEL



- Statistical design of experiments
- 29 experiments with different combinations of parameter values
- Experiment 1, 15, 29 “center points” with average parameter values
- CIS splitless time and further optimization with optimized TDU/CIS parameters

Parameter	Minimum	Maximum
TDU final temperature [°C]	250	300
TDU hold time [min]	3	10
TDU heating rate [°C/min]	80	250
Desorption flow [mL/min]	30	80
CIS trapping temperature [°C]	-30	-100


Method Optimization – TDU/CIS Parameters

GERSTEL




Trusted Answers


Benzo[g,h,i]perylene




Hexachlorobutadiene



Benzo[a]pyrene



Endosulfan



Experiment number

www.gerstel.com


Sample Prep Solution Method

Sample Preparation


2nd Extraction (17 h) at elevated temp.:
+ organic modifier



1st Extraction (5 h):
100 mL water sample
including ISTD mix



Thermal Desorption



Thermal Desorption Unit
(TDU, 90 – 300 °C)
Valveless Cooled Injection System
(CIS, -40 – 300 °C)
TDU Solvent Vent 2 min at 90°C

GC-MS/MS



7890/7010
HP-5ms UI 30m, 0.25 mm, 0.25 μm
60 °C, 1 min, 40 °C/min, 120 °C, 5
°C/min, 310 °C, 10 min
RTL chlorpyrifos-methyl
MRM

e.g. 100 pg/L in 100 mL sample → 10 pg per injection

Reduced Sample Volume!

GERSTEL



Standard water
analysis method
0.5 - 1 L



Twister method
0.1 L

GERSTEL



Agilent

Trusted Answers

How does the developed Method perform?



Limits of Quantification (LOQ)

Analyte	LOQ [ng/L]
Acenaphthene	1.0
Acenaphthylene	0.10
Aclonifen	0.56
Alachlor	0.42
Aldrin	0.066
Ametryn	0.069
Anthracene	0.061
Atrazine	0.18
Benz[a]anthracene	0.076
Benzo[a]pyrene	0.033
Benzo[b]fluoranthene	0.078
Benzo[g,h,i]perylene	0.049
Benzo[k]fluoranthene	0.081
Bifenox	0.47
Biphenyl	9.0
Bis(2-ethylhexyl)phthalate (DEHP)	134
Chlordane, cis	0.052
Chlordane, trans	0.026
Chlorfenvinphos	0.084
Chlorpyrifos-ethyl	0.024
Chrysene	0.027
Cybutryne (Irgarol 1051)	0.030
Cypermethrine (4 isomers)	0.12
p,p'-DDD	0.020
p,p'-DDE	0.017
o,p'-DDT	0.052
p,p'-DDT	0.067
Dibenz[a,h]anthracene	0.073
Dichlobenil	2.1
Dichlorvos	0.073
Dicofol	0.15
Dieldrin	0.034
Diflufenican	0.16
2,6-Di-tert.-butyl-4-methylphenol	5.9


Analyte	LOQ [ng/L]
alpha-Endosulfan	0.070
beta-Endosulfan	0.059
Endosulfan sulfate	0.052
Endrin	0.043
Endrin ketone	0.052
Ethofumesate	0.073
Fenitrothion	0.024
Fenpropimorph	0.13
Fluoranthene	1.0
Fluorene	0.45
Heptachlor	0.052
Heptachlor epoxide	0.052
Hexachlorobenzene	0.10
Hexachlorbutadiene	0.043
alpha-Hexachlorocyclohexane	0.052
beta-Hexachlorocyclohexane	0.13
gamma-Hexachlorocyclohexane	0.052
delta-Hexachlorocyclohexane	0.052
Indeno[1,2,3-cd]pyrene	0.044
Isodrin	0.16
Methoxychlor	0.083
Metolachlor	0.064
Naphthalene	5.0
Nonylphenol	8.8
Octylphenol	0.46
Oxadiazon	0.082
PBDE 28	0.018
PBDE 47	0.015
PBDE 99	0.050
PBDE 100	0.011
PBDE 153	0.032
PBDE 154	0.020
PBDE 183	0.13
PCB 77	0.041

Analyte	LOQ [ng/L]
PCB 81	0.039
PCB 105	0.043
PCB 114	0.036
PCB 118	0.012
PCB 123	0.037
PCB 126	0.050
PCB 156	0.046
PCB 157	0.047
PCB 167	0.044
PCB 169	0.054
PCB 189	0.054
Pendimethalin	0.094
Pentachlorobenzene	0.075
Pentachlorophenol	3.0
Phenanthrene	2.5
Picolinafen	0.26
Prometon	0.18
Prometryne	0.13
Propazine	0.057
Propiconazole	0.14
Propyzamide	0.35
Pyrene	0.45
Quinoxifen	0.087
Simazine	1.9
Terbutryne	0.1
Triallate	0.084
Tri-n-butyl phosphate	9.7
1,2,3-Trichlorobenzene	0.95
1,2,4-Trichlorobenzene	1.2
1,3,5-Trichlorobenzene	0.18
Triclosan	1.4
Trifluralin	0.19
Tris(2-chloroisopropyl)phosphate (TCPP)	29

LOQ determination according to DIN 32645 and under consideration of precision and trueness. If possible extracted from river water.


Uncertainty of measurements below 50 % for all EU-WFD compounds at or near the respective EQS as required by 2009/90/EC

Dichlorvos – River Water



Dichlorvos - 14 Levels, 14 Levels Used, 14 Points, 14 Points Used, 0 QCs

$y = 0.132458 * x - 0.002362$
 $R^2 = 0.99873142$
 Type:Linear, Origin:Ignore, Weight:1/x





0.073 ng/L in river water, LOQ 0.073 ng/L

7.3 pg on column (100 mL sample, assumption 100% recovery)

Required LOQ 0.18 ng/L (inland surface water)

Benzo[a]pyrene – Mineral Water





0.040 ng/L in mineral water, LOQ 0.033 ng/L

4 pg on column (100 mL sample, assumption 100% recovery)

Required LOQ 0.051 ng/L (inland surface water)

Cybutryne – Mineral Water





0.030 ng/L in mineral water, LOQ 0.030 ng/L

3 pg on column (100 mL sample, assumption 100% recovery)

Required LOQ 0.75 ng/L (inland surface water)

Cypermethrine – Ground Water



0.12 ng/L in ground water, LOQ **0.12 ng/L**

12 pg on column (100 mL sample, assumption 100% recovery)


Required LOQ **0.024 ng/L** (inland surface water)

Precision of Repeat Analyses [%]

GERSTEL




Spiked Mineral Water near LOQ (n=6)



Average RSD = 6.9 %

www.gerstel.com

Spiked Mineral Water near LOQ (n=6)



Average Deviation from 100 % = 7.2 %

Particle Load of Surface Water Samples*

- ▶ Data on suspended particulate matter (SPM) in characteristic French water samples, presented at CEN/TC 230 meeting – Paris, March 14-15, 2019
- ▶ 94% of water samples have a SPM load < 50 mg/L (with a 84% majority at SPM < 25 mg/L). Less than 1% of the samples had an SPM > 200 mg/L
- ▶ Effectiveness of LLE or SPE for extraction of analytes from SPM normally not proven/validated

* Personal communication Frank David, RIC


Extraction of Particle adsorbed Analytes



5, 10 or 15 mg
WEPAL SETOC 745
reference sediment




100 mL tap water
(50, 100 or 150 mg particles/L)
+ internal standard



Twister extraction under heating
TDU-GC-MS/MS
-> Analysis results

Extraction Parameters – “Sequential SBSE”

GERSTEL



- Expected analysis value
- No modifier, room temperature, 17 h
- 30% MeOH + 2 µL organic modifier, room temperature, 17 h
- 15% MeOH + 50 µL organic modifier, 65 °C, 17 h

Extraction of Particle adsorbed Analytes


10 mg Sediment in 100 mL Water (100 mg/L)

GERSTEL



Agilent


Trusted Answers



Extraction of Particle adsorbed Analytes

10 mg Sediment in 100 mL Water (100 mg/L)

GERSTEL



Externally spiked Ground Water Sample at 0.5 ng/L


GERSTEL



Trusted Answers



IWW WATER CENTRE



Analysis of real Water Samples

Concentrations in ng/L

GERSTEL



		Sewage Treatment Plant Effluent				Tap Water
Ametryne 3.4	Fluoranthene 0.73	Acenaphthylene 0.22	Dibenz[a,h]anthracene 0.20	4-Nonylphenols 199	Pyrene 2.5	Anthracene 0.069
Anthracene 0.28	alpha Hexachlorocyclohexane, 0.14	Alachlor 0.46	Dieldrin 0.084	4-Octylphenol 11	Terbutryne 6.1	trans-Chlordane 0.038
Atrazine 0.51	gamma Hexachlorocyclohexane 0.12	Anthracene 0.42	Diflufenican 0.39	PBDE 47 0.098	Triallate 0.16	alpha Hexachlorocyclohexane 0.065
Benz[a]anthracene 0.21	Indeno[1,2,3-cd]pyrene 0.39	Atrazine 1.6	2,6-Di-tert.-butyl-4-methylphenol 6.4	PBDE 99 0.23	Triclosan 9.9	beta Hexachlorocyclohexane 0.15
Benzo[a]pyrene 0.33	Metolachlor 8.2	Benz[a]anthracene 1.5	alpha Endosulfan 0.24	PBDE 100 0.046	Tri-n-butylphosphate 39	gamma Hexachlorocyclohexane 0.056
Benzo[b]fluoranthene 0.33	Nonylphenol 11	Benzo[a]pyrene 1.9	Ethofumesate 0.32	PBDE 154 0.031	Tris(2-chloroisopropyl) phosphate (TCPP) >300	4-Nonylphenols 5.9
Benzo[g,h,i]perylene 0.20	Propiconazol 0.57	Benzo[b]fluoranthene 2.1	Fluoranthene 2.3	PCB 118 0.021	42 compounds >LOQ	4-Octylphenol 2.3
Benzo[k]fluoranthene 0.20	Propyzamide 19	Benzo[g,h,i]perylene 2.0	alpha Hexachlorocyclohexane 0.071	Pendimethalin 0.36		7 compounds >LOQ
trans-Chlordane 0.035	Pyrene 0.53	Benzo[k]fluoranthene 1.1	beta Hexachlorocyclohexane, 0.31	Pentachlorobenzene 0.099		
Chrysene 0.40	Tebutryne 0.22	Chlorpyrifos-ethyl 0.23	gamma Hexachlorocyclohexane, 1.4	Propazine 0.11		
Cybutryne 1.1	Tris(2-chloroisopropyl) phosphate (TCPP) 52	Chrysene 1.9	Indeno[1,2,3-cd]pyrene 2.6	Propiconazol 6.1		
Diflufenican 2.6	23 compounds >LOQ	Cybutryne 0.28	Metolachlor 3.8	Propyzamide 2.7		

Conclusions

- ▶ Employing the **Twister** with **Agilent 7010 (HES)** is useful for EU-WFD guided water analysis
- ▶ **LOQ** of analytes mainly in the range of **two to three digit pg/L**
- ▶ Analysis method works for around **100 relevant compounds**, further can be added
- ▶ **Required LOQs** for inland surface water are achieved for all chosen analytes except heptachlor, heptachlor-epoxide, cypermethrine
- ▶ Extraction of **particle adsorbed compounds** is feasible and was confirmed for a reference sediment
- ▶ Ready to use GERSTEL **“Sample Prep Solution”** with defined hardware and method, manual, validation data etc.
- ▶ Many laboratories in Spain successfully use Twister for EU-WFD work, also labs in Germany (state lab in Augsburg, Bavaria and state lab in Neumünster, Schleswig Holstein) and other European countries

