



5TH SBSE INTERNATIONAL MEETING

23 & 24 SEPTEMBRE 2019 - NOVOTEL PARIS-SUD

SBSE 
Technical Meeting

ICE Concentration Linked with Extractive Stirrer (ICECLES)

SBSE “on the rocks” - David Benanou

Brian A. Logue

South Dakota State University

LARGE






Freeze Concentration (FC)

SOUTH DAKOTA
STATE UNIVERSITY

- Freeze concentration has been used for centuries, mainly in the food industry.
 - *Frozen concentrate orange juice*
 - *Ice beer*
 - *Ice wine*
- Sparingly evaluated for sample preparation.



Dia et al. (2011) *Environmental Engineering and Management Journal*




How does FC work?

SOUTH DAKOTA
STATE UNIVERSITY

- Freeze concentration works through freezing point depression.
- A higher concentration of solute in a small area of volume causes lower freezing points.
- Therefore, lower concentration areas of solution freeze first and increase the concentration of the resulting fluid.
- Under the right conditions, almost pure solvent can be frozen out.
- In our lab, we found that under the right conditions, the stir bar will also stay stirring on top of the ice...creating an obvious link to SBSE.

$$\Delta T_F = iK_F m$$




LARGE



ICECLES Apparatus

SOUTH DAKOTA
STATE UNIVERSITY



- Typical experimental...
 - Coolant at -3 to -7 °C.
 - Start at 1200 rpm; 300 rpm after 2 hr.
 - Aqueous sample.
 - Freezing over 2-3 hours.



LARGE



ICECLES Sample Prep

SOUTH DAKOTA
STATE UNIVERSITY



- ICECLES was performed using...
 - *1200 rpm stirring.*
 - *Coolant at -5 °C.*
 - *Methylene blue dye.*
 - *Freezing over 3 hours.*





LARGE



ICECLES Schematic

SOUTH DAKOTA
STATE UNIVERSITY



- White arrows designate the position of the stir bar.
- Stir bar is removed and typically...
 - *placed in a TDU tube for GC.*
 - *Back-extracted for LC.*



LARGE




ICECLES “Equilibrium”

SOUTH DAKOTA
STATE UNIVERSITY

SBSE f_{extr} :

$$f_{extr} = \frac{K_{SB}V_{SB}}{K_{SB}V_{SB} + V_{aq}}$$



ICECLES f_{extr} :

$$f_{extr} = \frac{K_{SB}V_{SB}}{K_{ice}V_{ice} + K_{SB}V_{SB} + V_{aq}} \rightarrow 0$$


- Need very small K_{ice} for large f_{extr} .
- Fortunately, K_{ice} is generally 10^{-2} to 10^{-4} .
- K_{ice} is a function of the rate of freezing, etc.





Fraction Extracted by ICECLES


SOUTH DAKOTA
STATE UNIVERSITY





Fraction Extracted by ICECLES


SOUTH DAKOTA
STATE UNIVERSITY





Signal Enhancement with ICECLES

SOUTH DAKOTA
STATE UNIVERSITY



Maslamani et al., (2016).




LARGE



Atrazine ($\log K_{ow} = 2.4$) Extraction

SOUTH DAKOTA
STATE UNIVERSITY





Extraction of Atrazine

SOUTH DAKOTA
STATE UNIVERSITY

Technique	Sample Size (mL)	Automated	Signal Relative to SBSE	Recovery (%)
SPE	1000	No	2.3	3.2
LLE	500	No	2.2	6.4
SBSE	10	Yes	1	14
SPME ¹	10	Yes	2.9×10^{-3}	ND ²
SPME ³	10	Yes	6.4×10^{-3}	ND ²
ICECLES	10	Yes	4.9	70

¹Back-extracted with 500 μ L

²ND = not detected


³Back-extracted with 55 μ L





Atrazine Extraction During Freezing


SOUTH DAKOTA
STATE UNIVERSITY





Analysis of NDPA in Drinking Water

SOUTH DAKOTA
STATE UNIVERSITY




- NDPA ($\log K_{ow} = 1.4$ est)
 - *Carcinogen.*
 - *EPA candidate contaminant.*
 - *EPA reporting limit = 7 ng/L.*
- Analysis of low ppt levels via a low-resolution GCMS





Analysis of NDPA in Drinking Water

SOUTH DAKOTA
STATE UNIVERSITY






- ICECLES allows detection on a low resolution GCMS (SIM).
- Quantification of NDPA at 2 ng/L was achievable in one drinking water source.





Pesticide Analysis in Drinking Water

SOUTH DAKOTA
STATE UNIVERSITY



- ICECLES detects more compounds at a greater sensitivity than SBSE and especially better than SPME.
 - *ICECLES, SBSE, SPME (55 μ L back-extract) were able to detect 53, 44, and 39 compounds, respectively at 10 ng/mL, and 32, 25, and 13 compounds, respectively, at 0.1 ng/mL.*

Skaggs, Alluhayb, and Logue (under review)





Pesticide Analysis In Drinking Water

SOUTH DAKOTA
STATE UNIVERSITY

Compound	Log K _{ow}	Concentration Detected (ng/mL)				Compound	Log K _{ow}	Concentration Detected (ng/mL)			
		ICECLES	SBSE	SPME ¹	SPME ²			ICECLES	SBSE	SPME ¹	SPME ²
Aldicarb	1.1	ND	ND	ND ³	ND	Fuberidazole	2.7	0.1	0.1	1	10
Aldicarb-sulfoxide	1.4	10	ND	ND	ND	Furathiocarb	4.7	0.1	1	1	10
Aldicarb-sulfone	1.4	10	ND	ND	ND	Iprovalicarb	3.2	0.1	0.1	0.1	1
Aminocarb	1.9	ND	ND	ND	ND	Isopropalin	1.4	10	10	10	ND
Amitraz	5.5	10	ND	ND	ND	Isocarbamid	2.0	0.1	0.1	1	10
Acibenzolar-S-Methyl	3.1	0.1	0.1	1	10	Mepronil	0.6	0.1	0.1	0.1	0.1
Bendiocarb	1.7	0.1	1	ND	10	Methiocarb	2.9	ND	ND	ND	ND
Butylate	4.2	0.1	1	10	ND	Methomyl	0.1	10	ND	ND	ND
Carbaryl	2.4	0.1	0.1	1	10	Metolcarb	1.7	0.1	1	10	ND
Carbendazim	1.5	10	ND	ND	ND	Napropamide	3.4	0.1	0.1	1	10
Carbofuran	2.3	1	1	10	ND	Naproanilide	4.4	0.1	0.1	0.1	1
Chlorantraniliprole	2.8	0.1	1	1	10	Oxamyl	-0.5	10	ND	ND	ND
Chlorodimeform	2.9	10	10	10	ND	Oxamyl oxime	0.2	ND	ND	ND	ND
Chlorbufam	3.6	0.1	0.1	1	ND	Oryzalin	3.7	1	ND	10	ND
Chlorpropham	3.8	0.1	0.1	1	10	Phenmedipham	2.7	1	10	ND	ND
Cymiazole	2.5	ND	ND	ND	ND	Pirimicarb	1.7	0.1	0.1	1	10
Cycloate	3.9	10	ND	ND	ND	Promecarb	3.1	0.1	0.1	0.1	1
Desmedipham	3.2	1	1	ND	ND	Propamocarb HCl	4.9	ND	ND	ND	ND
Diallate	3.3	1	1	1	10	Propanil	3.1	0.1	0.1	0.1	1
Dimepiperate	5.6	10	10	10	ND	Propham	2.6	0.1	1	10	ND
Dioxacarb	4.9	1	10	ND	ND	Thiabendazole	2.5	0.1	0.1	0.1	1
Diphenamid	2.4	0.1	0.1	0.1	1	Thiodicarb	1.6	0.1	1	1	10
EPTC	3.2	10	ND	ND	ND	Triallate	4.6	1	1	1	10
Ethiofencarb	2.0	ND	ND	ND	ND	Trichlamide	5.6	0.1	0.1	0.1	1
Etobenzanid	4.3	0.1	0.1	0.1	1	2,3,5-Trimethacarb	2.5	0.1	0.1	1	10
Fenfuram	2.6	0.1	0.1	1	10	3,4,5-Trimethacarb	2.6	0.1	0.1	0.1	1
Fenoxycarb	4.3	0.1	0.1	0.1	1	Vernolate	3.8	1	1	1	10
Fenthioicarb	3.3	0.1	0.1	0.1	1	XMC	2.3	0.1	0.1	1	10
Fenoxanil	3.5	0.1	0.1	0.1	1	Xylylcarb	2.1	1	10	ND	ND
Formetanate HCl	4.6	1	1	ND	ND	Zoxamide	3.8	0.1	0.1	1	10



Pesticide Analysis In Purified Water


SOUTH DAKOTA
STATE UNIVERSITY

Compound	Log K _{ow}	Concentration Detected (ng/mL)				Compound	Log K _{ow}	Concentration Detected (ng/mL)			
		ICECLES	SBSE	SPME ¹	SPME ²			ICECLES	SBSE	SPME ¹	SPME ²
Aldicarb	1.1	10	10	ND ³	ND	Fuberidazole	2.7	0.1	0.1	1	10
Aldicarb-sulfoxide	1.4	1	ND	ND	ND	Furathiocarb	4.7	0.1	0.1	0.1	1
Aldicarb-sulfone	1.4	1	ND	ND	ND	Iprovalicarb	3.2	0.1	0.1	0.1	1
Aminocarb	1.9	0.1	1	1	10	Isopropalin	1.4	10	10	10	10
Amitraz	5.5	10	10	ND	ND	Isocarbamid	2.0	0.1	10	10	ND
Acibenzolar-S-Methyl	3.1	0.1	0.1	1	10	Mepronil	0.6	0.1	0.1	0.1	1
Bendiocarb	1.7	0.1	1	1	10	Methiocarb	2.9	10	10	10	10
Butylate	4.2	0.1	0.1	0.1	1	Methomyl	0.1	0.1	ND	ND	ND
Carbaryl	2.4	0.1	0.1	1	10	Metolcarb	1.7	0.1	1	10	ND
Carbendazim	1.5	0.1	10	10	ND	Napropamide	3.4	0.1	0.1	0.1	0.1
Carbofuran	2.3	0.1	1	10	ND	Naproanilide	4.4	0.1	0.1	0.1	0.1
Chlorantraniliprole	2.8	0.1	1	10	10	Oxamyl	-0.5	10	ND	ND	ND
Chlorodimeform	2.9	0.1	0.1	1	ND	Oxamyl oxime	0.2	10	ND	ND	ND
Chlorbufam	3.6	0.1	0.1	1	ND	Oryzalin	3.7	0.1	1	10	ND
Chlorpropham	3.8	0.1	0.1	1	10	Phenmedipham	2.7	0.1	1	1	10
Cymiazole	2.5	0.1	0.1	0.1	1	Pirimicarb	1.7	0.1	0.1	1	10
Cycloate	3.9	0.1	0.1	0.1	1	Promecarb	3.1	0.1	0.1	0.1	1
Desmedipham	3.2	0.1	0.1	0.1	1	Propamocarb HCl	4.9	0.1	1	10	ND
Diallate	3.3	0.1	0.1	0.1	1	Propanil	3.1	0.1	0.1	1	10
Dimepiperate	5.6	0.1	0.1	0.1	0.1	Propham	2.6	0.1	1	10	ND
Dioxacarb	4.9	0.1	0.1	0.1	1	Thiabendazole	2.5	0.1	0.1	0.1	1
Diphenamid	2.4	0.1	0.1	0.1	1	Thiodicarb	1.6	0.1	1	10	ND
EPTC	3.2	0.1	0.1	1	10	Triallate	4.6	0.1	0.1	0.1	1
Ethiofencarb	2.0	10	10	10	ND	Trichlamide	5.6	0.1	0.1	0.1	1
Etobenzanid	4.3	0.1	0.1	0.1	1	2,3,5-Trimethacarb	2.5	0.1	1	1	10
Fenfuram	2.6	0.1	0.1	1	10	3,4,5-Trimethacarb	2.6	0.1	1	1	10
Fenoxycarb	4.3	0.1	0.1	0.1	0.1	Vernolate	3.8	0.1	0.1	0.1	1
Fenthioicarb	3.3	0.1	0.1	0.1	0.1	XMC	2.3	0.1	1	1	10
Fenoxanil	3.5	0.1	0.1	0.1	1	Xyllylcarb	2.1	10	10	ND	ND
Formetanate HCl	4.6	0.1	1	10	ND	Zoxamide	3.8	0.1	0.1	0.1	1



ICECLES of Green Tea

SOUTH DAKOTA
STATE UNIVERSITY



Alluhayb and Logue (2017).






ICECLES of Green Tea

SOUTH DAKOTA
STATE UNIVERSITY

Peak No./Category	Name	Odor	Log K _{ow}
<i>Alcohol</i>			
14	1-Pentanol	Fruity	1.33
15	2-Penten-1-ol, (Z)	Rubber	0.9*
84	Phenylethyl Alcohol	Rose	1.57
105	α -Terpineol (α , α -trimethyl 3-Cyclohexene-1-methanol)	Floral	3.28
<i>Heterocyclic</i>			
21	Methyl pyrazine	Nut	0.49
38	2,5-dimethyl pyrazine	Nut	1.03
<i>Aldehyde</i>			
22	Furfural	Caramel	0.83
37	(E,E)-2,4-Hexadienal	Citrus	1.37 ^c
111	5-(hydroxymethyl) furfural (5-(hydroxymethyl)-2-Furancarboxaldehyde)	Carmel	-0.09 ^c
<i>Ketone</i>			
82	Maltol	Caramel	0.02
91	Ketoisophorone (2,6,6-Trimethyl-2-cyclohexene-1,4-dione)	Floral	1 ^b
<i>Ester</i>			
95	Benzyl acetate	Fruit	1.96 ^b
66	γ -Undecalactone (5-heptyldihydro-2-(3H)-Furanone)	Fruit	0.7 ^b
<i>Phenol</i>			
136	Syringol (2,6-dimethoxy phenol)	Phenol	1.1 ^b
138	Eugenol	Clove	2.49

^blog K_{ow} values were calculated by using the difference between a log K_{ow} value of known compound and the query compound then estimated by an additive model with well-defined correction factors [26].

^clog K_{ow} values were calculated by using an atom/fragment contribution method via KOWWIN™ program [27].



Alluhayb and Logue (2017). Time (Min)



Potential of ICECLES

SOUTH DAKOTA
STATE UNIVERSITY

- Advantages:
 - *Can achieve extremely high concentration factors.*
 - *Well-suited for trace targeted analysis and comprehensive analysis.*
 - *Extends SBSE to compounds with $\log K_{ow} < 3$. This polarity range is difficult to analyze for most sample preparation techniques.*
 - *Well-suited for more thermally labile and volatile compounds.*
 - *Easier to back-extract compounds, making ICECLES more amenable for LC analysis.*
- Areas of improvement:
 - *Lengthy sample preparation times (depending on sample volume).*
 - *Simple implementation requires small sample volumes.*
 - *Precipitation of some compounds more likely under freezing temps.*
 - *Extra equipment, compared with SBSE, necessary to freeze the sample.*
 - *More difficult for low freezing point solvents.*

Addressed by automated
large volume system.





Acknowledgements

SOUTH DAKOTA
STATE UNIVERSITY


- LARGE group members.
- SACM and Quasim University.
- U.S. Joint Executive Office for Chem Bio Defense, Joint Program Management Protection W911SR-09-0059.
- Meeting organizers and participants.
- South Dakota winters, soda, and my laziness.





Questions???

SOUTH DAKOTA
STATE UNIVERSITY



brian.logue@sdstate.edu

<https://www.sdstate.edu/chemistry-biochemistry/brian-logue>

